



PRL Antibody

Product Code	CSB-PA12344A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01236
Immunogen	Recombinant Human Prolactin protein (29-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolactin (PRL), PRL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.