



PSMC6 Antibody, Biotin conjugated

Product Code	CSB-PA018896LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62333
Immunogen	Recombinant Human 26S proteasome regulatory subunit 10B protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome regulatory subunit 10B (26S proteasome AAA-ATPase subunit RPT4) (Proteasome 26S subunit ATPase 6) (Proteasome subunit p42), PSMC6, SUG2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMC6