



PSMD7 Antibody, Biotin conjugated

Product Code	CSB-PA018912LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51665
Immunogen	Recombinant Human 26S proteasome non-ATPase regulatory subunit 7 protein (1-324AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome non-ATPase regulatory subunit 7 (26S proteasome regulatory subunit RPN8) (26S proteasome regulatory subunit S12) (Mov34 protein homolog) (Proteasome subunit p40), PSMD7, MOV34L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMD7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.