



RASSF4 Antibody

Product Code	CSB-PA880956LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2L5
Immunogen	Recombinant Human Ras association domain-containing protein 4 protein (65-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras association domain-containing protein 4, RASSF4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RASSF4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.