





RBPMS Antibody

Product Code	CSB-PA04569A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93062
Immunogen	Recombinant Human RNA-binding protein with multiple splicing protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein with multiple splicing (RBP-MS) (Heart and RRM expressed sequence) (Hermes), RBPMS, HERMES
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RBPMS

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA04569A0Rb at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA04569A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.