




# Rgn Antibody

<b>Product Code</b>	CSB-PA019630LA01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q03336
<b>Immunogen</b>	Recombinant Rat Regucalcin protein (1-299AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat, Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Regucalcin (RC) (Gluconolactonase) (GNL) (EC 3.1.1.17) (Senescence marker protein 30) (SMP-30), Rgn, Smp30
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Others
<b>Target Names</b>	Rgn

## Image



Western blot  
 All lanes: Rgn antibody at 8µg/ml  
 Lane 1: A549 whole cell lysate  
 Lane 2: HCT116 whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 34 kDa  
 Observed band size: 34 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.