



# RNF138 Antibody, HRP conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA851971LB01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q8WVD3  |
| <b>Immunogen</b>           | Recombinant Human E3 ubiquitin-protein ligase RNF138 protein (2-245AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | HRP   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | E3 ubiquitin-protein ligase RNF138 (EC 2.3.2.27) (Nemo-like kinase-associated RING finger protein) (NLK-associated RING finger protein) (hNARF) (RING finger protein 138) (RING-type E3 ubiquitin transferase RNF138), RNF138, NARF |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling   |
| <b>Target Names</b>        | RNF138  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |