



RPGR Antibody, HRP conjugated

Product Code	CSB-PA838814LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92834
Immunogen	Recombinant Human X-linked retinitis pigmentosa GTPase regulator protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	X-linked retinitis pigmentosa GTPase regulator, RPGR, RP3 XLRP3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RPGR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.