




RPS24 Antibody

Product Code	CSB-PA020402LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62847
Immunogen	Recombinant Human 40S ribosomal protein S24 protein (2-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S24 (Small ribosomal subunit protein eS24), RPS24
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPS24

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA020402LA01HU at dilution of 1:100