



RPS28 Antibody, HRP conjugated

Product Code	CSB-PA020430LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62857
Immunogen	Recombinant Human 40S ribosomal protein S28 protein (1-69AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S28 (Small ribosomal subunit protein eS28), RPS28
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS28
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.