



S100A11 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA020624HB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P31949 |
| Immunogen | Recombinant Human Protein S100-A11 protein (2-105AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein S100-A11 (Calgizzarin) (Metastatic lymph node gene 70 protein) (MLN 70) (Protein S100-C) (S100 calcium-binding protein A11) [Cleaved into: Protein S100-A11, N-terminally processed], S100A11, MLN70 S100C |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | S100A11 |