



SEC61G Antibody, FITC conjugated

Product Code	CSB-PA03509C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60059
Immunogen	Recombinant Human Protein transport protein Sec61 subunit gamma protein (1-32AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein transport protein Sec61 subunit gamma, SEC61G
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SEC61G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.