



SGTB Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA846614LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96EQ0 |
| Immunogen | Recombinant Human Small glutamine-rich tetratricopeptide repeat-containing protein beta protein (14-304AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Small glutamine-rich tetratricopeptide repeat-containing protein beta (Beta-SGT) (Small glutamine-rich protein with tetratricopeptide repeats 2), SGTB, SGT2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | SGTB |