



SH3BGRL Antibody, Biotin conjugated

Product Code	CSB-PA021220LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75368
Immunogen	Recombinant Human SH3 domain-binding glutamic acid-rich-like protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-binding glutamic acid-rich-like protein, SH3BGRL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SH3BGRL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.