



SLC25A25 Antibody, Biotin conjugated

Product Code	CSB-PA760850LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6KCM7
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-2 protein (1-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein SCaMC-2 (Mitochondrial ATP-Mg/Pi carrier protein 3) (Mitochondrial Ca(2+)-dependent solute carrier protein 3) (Small calcium-binding mitochondrial carrier protein 2) (Solute carrier family 25 member 25), SLC25A25, APC3 KIAA1896 MCSC3 SCAMC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A25