



# SLC25A25 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA760850LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6KCM7
<b>Immunogen</b>	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-2 protein (1-189AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Calcium-binding mitochondrial carrier protein SCaMC-2 (Mitochondrial ATP-Mg/Pi carrier protein 3) (Mitochondrial Ca(2+)-dependent solute carrier protein 3) (Small calcium-binding mitochondrial carrier protein 2) (Solute carrier family 25 member 25), SLC25A25, APC3 KIAA1896 MCSC3 SCAMC2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SLC25A25
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.