



SLC25A40 Antibody, Biotin conjugated

Product Code	CSB-PA840561LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TBP6
Immunogen	Recombinant Human Solute carrier family 25 member 40 protein (40-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 25 member 40 (Mitochondrial carrier family protein), SLC25A40, MCFP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC25A40