



SLC39A11 Antibody, Biotin conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA854023LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8N1S5 |
| Immunogen | Recombinant Human Zinc transporter ZIP11 protein (93-193AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc transporter ZIP11 (Solute carrier family 39 member 11) (Zrt- and Irt-like protein 11) (ZIP-11), SLC39A11, C17orf26 ZIP11 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | SLC39A11 |