



SLC4A8 Antibody, Biotin conjugated

Product Code	CSB-PA649721LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2Y0W8
Immunogen	Recombinant Human Electroneutral sodium bicarbonate exchanger 1 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Electroneutral sodium bicarbonate exchanger 1 (Electroneutral Na(+)-driven Cl-HCO ₃ exchanger) (Solute carrier family 4 member 8) (k-NBC3), SLC4A8, KIAA0739 NBC NBC3 NDCBE1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC4A8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.