




SPATS2L Antibody

Product Code	CSB-PA868324HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUQ6
Immunogen	Recombinant Human SPATS2-like protein (2-488AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SPATS2-like protein (DNA polymerase-transactivated protein 6) (Stress granule and nucleolar protein) (SGNP), SPATS2L, DNAPTP6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPATS2L

Image



Western blot

All lanes: SPATS2L antibody at 8µg/ml

Lane 1: CEM whole cell lysate

Lane 2: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 62, 55, 66, 63 kDa

Observed band size: 62 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.