



STEAP3 Antibody, Biotin conjugated

Product Code	CSB-PA717655LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q658P3
Immunogen	Recombinant Human Metalloreductase STEAP3 protein (1-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metalloreductase STEAP3 (EC 1.16.1.-) (Dudulin-2) (Six-transmembrane epithelial antigen of prostate 3) (Tumor suppressor-activated pathway protein 6) (hTSAP6) (pHyde) (hpHyde), STEAP3, TSAP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STEAP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.