



SV2B Antibody, HRP conjugated

Product Code	CSB-PA022979LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L1I2
Immunogen	Recombinant Human Synaptic vesicle glycoprotein 2B protein (1-108AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptic vesicle glycoprotein 2B, SV2B, KIAA0735
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SV2B