




TBC1D22A Antibody

Product Code	CSB-PA848818LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUA7
Immunogen	Recombinant Human TBC1 domain family member 22A protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 22A, TBC1D22A, C22orf4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D22A

Image



Western blot

All lanes: TBC1D22A antibody at 18µg/ml

Lane 1: Colo320 whole cell lysate

Lane 2: HepG2 whole cell lysate

Lane 3: Rat skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 60, 51, 37, 54 kDa

Observed band size: 60 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.