



TBCB Antibody, FITC conjugated

Product Code	CSB-PA03169C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99426
Immunogen	Recombinant Human Tubulin-folding cofactor B protein (1-244AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin-folding cofactor B (Cytoskeleton-associated protein 1) (Cytoskeleton-associated protein CKAPI) (Tubulin-specific chaperone B), TBCB, CG22 CKAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBCB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.