



TBX15 Antibody, Biotin conjugated

Product Code	CSB-PA839414HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96SF7
Immunogen	Recombinant Human T-box transcription factor TBX15 protein (1-496AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-box transcription factor TBX15 (T-box protein 15) (T-box transcription factor TBX14) (T-box protein 14), TBX15, TBX14
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TBX15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.