



TFEC Antibody

Product Code	CSB-PA023430LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14948
Immunogen	Recombinant Human Transcription factor EC protein (1-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor EC (TFE-C) (Class E basic helix-loop-helix protein 34) (bHLHe34) (Transcription factor EC-like) (hTFEC-L), TFEC, BHLHE34 TC FEC TFEC L
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TFEC

Image



Western blot

All lanes: TFEC antibody at 3µg/ml + MDA-MB-231 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 39, 36, 23, 32 kDa

Observed band size: 39 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.