



TGFB3 Antibody, Biotin conjugated

Product Code	CSB-PA023449LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10600
Immunogen	Recombinant Human Transforming growth factor beta-3 proprotein protein (24-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-3 proprotein (TGF-beta-3) [Cleaved into: Latency-associated peptide (LAP)] (Transforming growth factor beta-3), TGFB3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TGFB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.