



TM4SF1 Antibody, Biotin conjugated

Product Code	CSB-PA023615LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30408
Immunogen	Recombinant Human Transmembrane 4 L6 family member 1 protein (115-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane 4 L6 family member 1 (Membrane component chromosome 3 surface marker 1) (Tumor-associated antigen L6), TM4SF1, M3S1 TAAL6
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	TM4SF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.