



TMED3 Antibody, Biotin conjugated

Product Code	CSB-PA897107LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3Q3
Immunogen	Recombinant Human Transmembrane emp24 domain-containing protein 3 protein (24-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane emp24 domain-containing protein 3 (Membrane protein p24B) (p24 family protein gamma-4) (p24gamma4) (p26), TMED3, C15orf22
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMED3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.