



TMEM115 Antibody, HRP conjugated

Product Code	CSB-PA023683LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12893
Immunogen	Recombinant Human Transmembrane protein 115 protein (187-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 115 (Placental protein 6) (Protein PL6), TMEM115, PL6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMEM115
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.