



TMEM59L Antibody, Biotin conjugated

Product Code	CSB-PA023861LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK28
Immunogen	Recombinant Human Transmembrane protein 59-like protein (27-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 59-like (Brain-specific membrane-anchored protein), TMEM59L, BSMAP C19orf4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TMEM59L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.