



TMEM61 Antibody, FITC conjugated

Product Code	CSB-PA023863LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0U2
Immunogen	Recombinant Human Transmembrane protein 61 protein (90-210AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 61, TMEM61
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM61