



TRIML1 Antibody, Biotin conjugated

Product Code	CSB-PA822741LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9V2
Immunogen	Recombinant Human Probable E3 ubiquitin-protein ligase TRIML1 protein (277-468AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable E3 ubiquitin-protein ligase TRIML1 (EC 2.3.2.27) (RING finger protein 209) (RING-type E3 ubiquitin transferase TRIML1) (Tripartite motif family-like protein 1), TRIML1, RNF209
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRIML1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.