





TSPAN14 Antibody

Product Code	CSB-PA847712LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NG11
Immunogen	Recombinant Human Tetraspanin-14 protein (114-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-14 (Tspan-14) (DC-TM4F2) (Transmembrane 4 superfamily member 14), TSPAN14, TM4SF14
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TSPAN14


Image



Immunofluorescent analysis of U251 cells using CSB-PA847712LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA847712LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA847712LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.