



TSPAN15 Antibody, FITC conjugated

Product Code	CSB-PA025152LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95858
Immunogen	Recombinant Human Tetraspanin-15 protein (115-235AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-15 (Tspan-15) (Tetraspan NET-7) (Transmembrane 4 superfamily member 15), TSPAN15, NET7 TM4SF15
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TSPAN15