





TSPAN5 Antibody

Product Code	CSB-PA025163LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62079
Immunogen	Recombinant Human Tetraspanin-5 protein (114-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-5 (Tspan-5) (Tetraspan NET-4) (Transmembrane 4 superfamily member 9), TSPAN5, TM4SF9
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TSPAN5

Image



Western blot
 All lanes: TSPAN5 antibody at 1.2µg/ml + Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 31 kDa
 Observed band size: 31 kDa



Immunofluorescent analysis of MCF-7 cells using CSB-PA025163LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.