



TTLL5 Antibody, Biotin conjugated

Product Code	CSB-PA715737LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6EMB2
Immunogen	Recombinant Human Tubulin polyglutamylase TTLL5 protein (1-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin polyglutamylase TTLL5 (EC 6.-.-.-) (SRC1 and TIF2-associated modulatory protein) (Tubulin--tyrosine ligase-like protein 5), TTLL5, KIAA0998 STAMP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTLL5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.