



UBA7 Antibody, Biotin conjugated

Product Code	CSB-PA025421LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41226
Immunogen	Recombinant Human Ubiquitin-like modifier-activating enzyme 7 protein (1-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like modifier-activating enzyme 7 (Ubiquitin-activating enzyme 7) (D8) (Ubiquitin-activating enzyme E1 homolog), UBA7, UBE1L UBE2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBA7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.