



UBLCP1 Antibody, HRP conjugated

Product Code	CSB-PA851974LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVY7
Immunogen	Recombinant Human Ubiquitin-like domain-containing CTD phosphatase 1 protein (2-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like domain-containing CTD phosphatase 1 (EC 3.1.3.16) (Nuclear proteasome inhibitor UBLCP1), UBLCP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBLCP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.