



UHRF1BP1L Antibody

Product Code	CSB-PA025606LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0JNW5
Immunogen	Recombinant Human UHRF1-binding protein 1-like protein (251-522AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UHRF1-binding protein 1-like, UHRF1BP1L, KIAA0701
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UHRF1BP1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.