





UNC50 Antibody

Product Code	CSB-PA687496LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53HI1
Immunogen	Recombinant Human Protein unc-50 homolog protein (15-82AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein unc-50 homolog (Periodontal ligament-specific protein 22) (PDLs22) (Protein GMH1 homolog) (hGMH1) (Uncoordinated-like protein), UNC50, UNCL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UNC50

Image



Western blot
 All lanes: UNC50 antibody at 0.2µg/ml + 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 31 kDa
 Observed band size: 31 kDa



Immunofluorescent analysis of A549 cells using CSB-PA687496LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.