



UXS1 Antibody, Biotin conjugated

Product Code	CSB-PA822754LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBZ7
Immunogen	Recombinant Human UDP-glucuronic acid decarboxylase 1 protein (201-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glucuronic acid decarboxylase 1 (EC 4.1.1.35) (UDP-glucuronate decarboxylase 1) (UGD) (UXS-1), UXS1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UXS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.