



VIPAS39 Antibody, HRP conjugated

Product Code	CSB-PA888005LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9C1
Immunogen	Recombinant Human Spermatogenesis-defective protein 39 homolog protein (1-493AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatogenesis-defective protein 39 homolog (hSPE-39) (VPS33B-interacting protein in apical-basolateral polarity regulator) (VPS33B-interacting protein in polarity and apical restriction), VIPAS39, C14orf133 SPE39 VIPAR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VIPAS39