



VIPAS39 Antibody, HRP conjugated

Product Code	CSB-PA888005LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9C1
Immunogen	Recombinant Human Spermatogenesis-defective protein 39 homolog protein (1-493AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatogenesis-defective protein 39 homolog (hSPE-39) (VPS33B-interacting protein in apical-basolateral polarity regulator) (VPS33B-interacting protein in polarity and apical restriction), VIPAS39, C14orf133 SPE39 VIPAR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VIPAS39
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.