



VPS26B Antibody, Biotin conjugated

Product Code	CSB-PA676824LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4G0F5
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 26B protein (1-336AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 26B (Vesicle protein sorting 26B), VPS26B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VPS26B