



VTA1 Antibody, Biotin conjugated

Product Code	CSB-PA882062LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP79
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein VTA1 homolog protein (23-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein VTA1 homolog (Dopamine-responsive gene 1 protein) (DRG-1) (LYST-interacting protein 5) (LIP5) (SKD1-binding protein 1) (SBP1), VTA1, C6orf55
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	VTA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.