



VWA3B Antibody, Biotin conjugated

Product Code	CSB-PA681377LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q502W6
Immunogen	Recombinant Human von Willebrand factor A domain-containing protein 3B protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	von Willebrand factor A domain-containing protein 3B (VWA domain-containing protein 3B), VWA3B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VWA3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.