



WDR88 Antibody, FITC conjugated

Product Code	CSB-PA744395LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZMY6
Immunogen	Recombinant Human WD repeat-containing protein 88 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 88 (PQQ repeat and WD repeat-containing protein), WDR88, PQWD
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDR88
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.