



ZBTB32 Antibody, Biotin conjugated

Product Code	CSB-PA897538LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2Y4
Immunogen	Recombinant Human Zinc finger and BTB domain-containing protein 32 protein (1-294AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and BTB domain-containing protein 32 (FANCC-interacting protein) (Fanconi anemia zinc finger protein) (Testis zinc finger protein) (Zinc finger protein 538), ZBTB32, FAZF TZFP ZNF538
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZBTB32
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.