




ZC3H7A Antibody

Product Code	CSB-PA818245LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWR0
Immunogen	Recombinant Human Zinc finger CCCH domain-containing protein 7A protein (1-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger CCCH domain-containing protein 7A, ZC3H7A, ZC3H7 ZC3HDC7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZC3H7A

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA818245LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.